

Product datasheet for **HP226188**

GCNF (NR6A1) Human qPCR Primer Pair (NM_033334)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	2649
Forward Sequence:	GGCTGATGTCACCTGCCAAGTAC
Reverse Sequence:	GCATACTCCTCGTTGCTGACCT
ACCN:	<u>NM_033334</u> , <u>NM_033334.1</u> , <u>NM_033334.2</u> , <u>NM_033334.3</u> , <u>BC030600</u> , <u>BC030600.1</u> , <u>BC027988</u> , <u>BC042395</u> , <u>BC042406</u> , <u>NM_033334.4</u>
UniProt ID:	<u>Q15406</u>
Synonyms:	CT150; GCNF; GCNF1; hGCNF; hRTR; NR61; RTR
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60 °C for 1 min; Stage 4: Melting curve: 95 °C for 15 sec, 60 °C for 15 sec, 95 °C for 15 sec.
Storage:	Store at -20 °C.
Stability:	The primer mix is stable for one year from date of shipping.


[View online »](#)