

Product datasheet for **HP226233**

IL22 RA2 (IL22RA2) Human qPCR Primer Pair (NM_181309)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	116379
Forward Sequence:	AGGGCACTTACTGGCAACAGCA
Reverse Sequence:	TGTCTGAGGTTTCACTGGTAAGG
ACCN:	<u>NM_181309</u> , <u>NM_181309.1</u> , <u>BC125167</u> , <u>BC125168</u> , <u>NM_181309.2</u>
UniProt ID:	<u>Q969J5</u>
Synonyms:	CRF2-10; CRF2-S1; CRF2X; IL-22BP; IL-22R-alpha-2; IL-22RA2; ZCYTOR16
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)