

Product datasheet for **HP226241**

TAF A3 Human qPCR Primer Pair (NM_182759)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	284467
Forward Sequence:	CCAGACGGTGAAATGCTCCTGT
Reverse Sequence:	AGGCAGGGCTCCATCTGACAC
ACCN:	NM_182759 , NM_182759.1 , NM_182759.2 , BC117479 , BC112268 , BC143949 , BC144042 , BC144043 , BQ638101
UniProt ID:	Q7Z5A8
Synonyms:	FAM19A3; TAF A-3
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)