

Product datasheet for **HP226258**

PORCN Human qPCR Primer Pair (NM_203475)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	64840
Forward Sequence:	TGGATGACACCACAGAGGAGCA
Reverse Sequence:	CCTCAGCCTATGAGACGGTAGA
ACCN:	<u>NM_203475</u> , <u>NM_203475.1</u> , <u>NM_203475.2</u> , <u>BC019080</u> , <u>BC019080.2</u> , <u>BM543485</u> , <u>BM545893</u> , <u>NM_203475.3</u>
UniProt ID:	<u>Q9H237</u>
Synonyms:	DHOF; FODH; MG61; PORC; PPN
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)