

Product datasheet for **HP226270**

GABA A Receptor beta 2 (GABRB2) Human qPCR Primer Pair (NM_021911)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	2561
Forward Sequence:	GTGCCAACAAATGAGAAGATGCGC
Reverse Sequence:	TTCCAGTAGGGTCCCACAAGGA
ACCN:	<u>NM_021911</u> , <u>NM_021911.1</u> , <u>NM_021911.2</u> , <u>BC099705</u> , <u>BC036592</u> , <u>BC099719</u> , <u>BC105639</u>
UniProt ID:	<u>P47870</u>
Synonyms:	DEE92; ICEE2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)