

Product datasheet for **HP226271**

Caspase 1 (CASP1) Human qPCR Primer Pair (NM_033292)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	834
Forward Sequence:	GCTGAGGTTGACATCACAGGCA
Reverse Sequence:	TGCTGTCAGAGGTCTTGCTC
ACCN:	<u>NM_033292</u> , <u>NM_033292.1</u> , <u>NM_033292.2</u> , <u>NM_033292.3</u> , <u>BC062327</u> , <u>BC062327.1</u> , <u>BC041689</u> , <u>NM_033292.4</u>
UniProt ID:	<u>P29466</u>
Synonyms:	ICE; IL1BC; P45
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)