

Product datasheet for **HP226275**

Rho guanine exchange factor 15 (ARHGEF15) Human qPCR Primer Pair (NM_173728)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	22899
Forward Sequence:	CACTCTTCTCCAATGTGCAGCG
Reverse Sequence:	CCGCACATAATCCACATACACCG
ACCN:	<u>NM_173728</u> , <u>NM_173728.1</u> , <u>NM_173728.2</u> , <u>NM_173728.3</u> , <u>BC036749</u> , <u>BC036749.1</u> , <u>BQ004849</u> , <u>BQ877184</u> , <u>NM_173728.4</u>
UniProt ID:	<u>O94989</u>
Synonyms:	ARGEF15; E5; Ephexin5; Vsm-RhoGEF
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)