

## Product datasheet for **HP226290**

### CCNK Human qPCR Primer Pair (NM\_003858)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	8812
Forward Sequence:	CACTATGATACCCTGGCAACTGG
Reverse Sequence:	CAGAAAGAGGCAACAGGCTCCT
ACCN:	<a href="#">NM_003858</a> , <a href="#">NM_003858.2</a> , <a href="#">BC015935</a> , <a href="#">BG503760</a> , <a href="#">BM505592</a> , <a href="#">BQ432222</a> , <a href="#">BT006950</a> , <a href="#">BX247958</a>
UniProt ID:	<a href="#">O75909</a>
Synonyms:	CPR4; cyclin K; MGC9113
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)