

Product datasheet for **HP226291**

RGR Human qPCR Primer Pair (NM_002921)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	5995
Forward Sequence:	GACAGTGGGATCAGCCTGAATG
Reverse Sequence:	AACGCTGTCACAAAGCCCTGGA
ACCN:	BC008094 , NM_002921 , NM_002921.1 , NM_002921.2 , NM_002921.3 , BC008094.1 , BC011349 , BC011349.1 , BC027987 , BC042536 , BG912392 , BX641033
UniProt ID:	P47804
Synonyms:	RP44
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)