

Product datasheet for **HP226394**

SLC35C2 Human qPCR Primer Pair (NM_173073)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	51006
Forward Sequence:	TCGTTTCATCGGTGGCATTTCGCT
Reverse Sequence:	TGCAGGTGGAACATGGTGTCTGA
ACCN:	<u>NM_173073</u> , <u>NM_173073.1</u> , <u>NM_173073.2</u> , <u>NM_173073.3</u> , <u>BC014191</u> , <u>BC025277</u> , <u>BC025277.1</u> , <u>BC021138</u> , <u>BI831689</u> , <u>BM666377</u> , <u>BM846744</u> , <u>NM_173073.4</u>
UniProt ID:	<u>Q9NQQ7</u>
Synonyms:	BA394O2.1; C20orf5; CGI-15; OVCOV1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.


[View online »](#)