

## Product datasheet for **HP226411**

### C 4 Methylsterol Oxidase (MSMO1) Human qPCR Primer Pair (NM\_006745)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	6307
Forward Sequence:	GCTGCCTTTGATTTGTGGAACCT
Reverse Sequence:	CTGCACAACCAAAGCATCTTGCC
ACCN:	<a href="#">NM_006745</a> , <a href="#">NM_006745.1</a> , <a href="#">NM_006745.2</a> , <a href="#">NM_006745.3</a> , <a href="#">NM_006745.4</a> , <a href="#">BC107879</a> , <a href="#">BC107879.1</a> , <a href="#">BC010653</a> , <a href="#">BU676028</a> , <a href="#">BU940384</a> , <a href="#">BX441001</a> , <a href="#">NM_006745.5</a>
UniProt ID:	<a href="#">Q15800</a>
Synonyms:	DESP4; ERG25; MCCPD; SC4MOL
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)