

Product datasheet for **HP226419**

Retinoid X Receptor gamma (RXRG) Human qPCR Primer Pair (NM_006917)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	6258
Forward Sequence:	GACTGCCTCATTGACAAGCGTC
Reverse Sequence:	GACCACTGGTAGCACATTCTGC
ACCN:	NM_006917 , NM_006917.1 , NM_006917.2 , NM_006917.3 , NM_006917.4 , BC012063 , BC012063.1 , BE349162 , BI768747 , BM684530 , BM728015 , BX118484 , NM_006917.5
UniProt ID:	P48443
Synonyms:	NR2B3; RXRC
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)