

Product datasheet for **HP226461**

PHD4 (P4HTM) Human qPCR Primer Pair (NM_177938)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	54681
Forward Sequence:	TGGTGGATGACTCCAGAGAGCA
Reverse Sequence:	ACTCTGCCTTGTGGCTCCTCAT
ACCN:	<u>NM_177938</u> , <u>NM_177938.1</u> , <u>NM_177938.2</u> , <u>BC011710</u> , <u>BC047566</u> , <u>BC047566.1</u> , <u>BC000580</u> , <u>BC018731</u> , <u>BC060321</u>
UniProt ID:	<u>Q9NXG6</u>
Synonyms:	EGLN4; HIDEA; HIFPH4; P4H-TM; PH-4; PH4; PHD4
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)