

Product datasheet for **HP226462**

LYPD1 Human qPCR Primer Pair (NM_144586)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	116372
Forward Sequence:	GCTGCAAATCCAGTGCTACCAG
Reverse Sequence:	TGATGCACAGGACTTGCGGTAC
ACCN:	<u>NM_144586</u> , <u>NM_144586.1</u> , <u>NM_144586.2</u> , <u>NM_144586.3</u> , <u>NM_144586.4</u> , <u>NM_144586.5</u> , <u>NM_144586.6</u> , <u>BC017318</u> , <u>BC017318.2</u> , <u>BC011449</u> , <u>BC040046</u> , <u>BE386060</u> , <u>BG057775</u> , <u>BP276862</u> , <u>BY800190</u> , <u>NM_144586.7</u>
UniProt ID:	<u>Q8N2G4</u>
Synonyms:	LYPDC1; PHTS
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)