

Product datasheet for **HP226472**

PEX16 Human qPCR Primer Pair (NM_057174)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	9409
Forward Sequence:	TCTTGGCTGGTGTGTGGACGT
Reverse Sequence:	AAAGGAGAGCGCAGCAGGTAGT
Accession No:	<u>NM_004813</u> , <u>NM_057174</u> , <u>BC004356</u> , <u>BC000467</u> , <u>BM746616</u> , <u>NM_057174.3</u>
Synonyms:	PBD8A; PBD8B
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	The primer mix is stable for one year from date of shipping. Store at -20°C.



[View online »](#)