

Product datasheet for **HP226482**

PRMT2 Human qPCR Primer Pair (NM_206962)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	3275
Forward Sequence:	GCAGTTGGACATGAGAACCGTG
Reverse Sequence:	AGGCTCTGGAAGTGGACGCTAA
ACCN:	NM_206962 , NM_206962.1 , NM_206962.2 , NM_206962.3 , BC000727 , BC038765 , BC100026 , BF683143 , BM976817 , BQ267811 , BQ421001 , BY799824 , NM_206962.4
UniProt ID:	P55345
Synonyms:	HRMT1L1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)