

Product datasheet for **HP226535**

PRR4 Human qPCR Primer Pair (NM_007244)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	11272
Forward Sequence:	CCCGATTTCTTCTGTGTCAGCCT
Reverse Sequence:	TGGCTTTCTGAAGGAAGTTATCTTC
ACCN:	<u>NM_007244</u> , <u>NM_007244.1</u> , <u>NM_007244.2</u> , <u>BC058035</u> , <u>BC058035.1</u> , <u>BX537839</u>
UniProt ID:	<u>Q16378</u>
Synonyms:	LPRP; PROL4
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)