

Product datasheet for **HP226554**

GALK2 Human qPCR Primer Pair (NM_002044)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	2585
Forward Sequence:	GCTTTGGTCTGTTGTGCTGGCT
Reverse Sequence:	TGTAACGCTCACTCTGGCACAG
ACCN:	NM_002044 , NM_002044.1 , NM_002044.2 , NM_002044.3 , BC005141 , BC005141.2 , BI766538 , BQ921903 , BT006901 , BX096260 , BX455251
UniProt ID:	Q01415
Synonyms:	GK2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)