

Product datasheet for **HP226577**

DPH3 Human qPCR Primer Pair (NM_206831)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	285381
Forward Sequence:	TCTATCCCTGCCCATGTGGAGA
Reverse Sequence:	CTGGGACTGTTTCTCCACACAC
ACCN:	BC010181 , NM_206831 , NM_206831.1 , NM_206831.2 , BC010181.2 , BC073931 , BC094863 , NM_206831.3
UniProt ID:	Q96FX2
Synonyms:	DELGIP; DELGIP1; DESR1; DPH3A; KTI11; ZCSL2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)