

## Product datasheet for **HP226639**

### ethanolamine kinase (ETNK1) Human qPCR Primer Pair (NM\_018638)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	55500
Forward Sequence:	CCACAGGATTTGCAGATGAAGAC
Reverse Sequence:	TGGCAAAGCACAAACAGGTGAGC
ACCN:	<a href="#">BC006111</a> , <a href="#">NM_018638</a> , <a href="#">NM_018638.1</a> , <a href="#">NM_018638.2</a> , <a href="#">NM_018638.3</a> , <a href="#">NM_018638.4</a> , <a href="#">BC006111.1</a> , <a href="#">BC066907</a> , <a href="#">BC066907.1</a> , <a href="#">BC037324</a> , <a href="#">BC045794</a> , <a href="#">BT007084</a>
UniProt ID:	<a href="#">Q9HBU6</a>
Synonyms:	EKI; EK1; EK11; Nbla10396
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)