

Product datasheet for **HP226653**

XBP1 Human qPCR Primer Pair (NM_005080)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	7494
Forward Sequence:	CTGCCAGAGATCGAAAGAAGGC
Reverse Sequence:	CTCCTGGTTCTCAACTACAAGGC
ACCN:	<u>NM_005080</u> , <u>NM_005080.1</u> , <u>NM_005080.2</u> , <u>NM_005080.3</u> , <u>BC012841</u> , <u>BC012841.1</u> , <u>BC000938</u> , <u>BC015709</u> , <u>BI545390</u>
UniProt ID:	<u>P17861</u>
Synonyms:	TREB-5; TREB5; XBP-1; XBP2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)