

Product datasheet for **HP226678**

TORC1 (CRTC1) Human qPCR Primer Pair (NM_015321)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	23373
Forward Sequence:	CGCGTACTATGAGCAGCAGATG
Reverse Sequence:	CTGCGAGTAGTTGAGTGTGGAG
ACCN:	<u>NM_015321</u> , <u>NM_015321.1</u> , <u>NM_015321.2</u> , <u>BC028050</u> , <u>BC028050.1</u> , <u>BC017075</u> , <u>BC023614</u> , <u>BM680095</u>
UniProt ID:	<u>Q6UUV9</u>
Synonyms:	Mam-2; MAML2; MECT1; TORC-1; TORC1; WAMTP1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)