

Product datasheet for **HP226703**

ILRUN Human qPCR Primer Pair (NM_024294)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	64771
Forward Sequence:	GTGCTCATCTCCGAGTTCCAGA
Reverse Sequence:	GGCACACTGATGTTTGGGCTCT
ACCN:	NM_024294 , NM_024294.2 , NM_024294.3 , BC075810 , BC075810.1 , BC002328 , BC010143 , BC010184 , BC013241 , BC071885 , BC107732 , BE898102 , NM_024294.4
UniProt ID:	Q9H6K1
Synonyms:	C6orf106; dJ391O22.4; FP852
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)