

Product datasheet for **HP226747**

TCL6 Human qPCR Primer Pair (NM_014418)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	27004
Forward Sequence:	TCGGACGAAGTGTTCAGACC
Reverse Sequence:	TGCAGAGGCTCTCAAGCTGTAG
ACCN:	<u>NM_014418</u> , <u>NM_014418.1</u> , <u>NM_014418.2</u> , <u>NM_014418.3</u> , <u>BC041075</u> , <u>BC041075.1</u> , <u>BX161462</u>
Synonyms:	OTTHUMP00000195767; T-cell leukemia/lymphoma 6; T-cell leukemia/lymphoma 6 ORF72; T-cell leukemia/lymphoma 6 ORF141; TCL1-neighboring gene 1; TCL1-neighboring gene 2; TCL6e1; TCL6f1; TCL6f2; TNG1; TNG2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)