

## Product datasheet for **HP226771**

### Macrophage Scavenger Receptor I (MSR1) Human qPCR Primer Pair (NM\_138715)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	4481
Forward Sequence:	TGCACAAGGCAGCTCACTTTGG
Reverse Sequence:	GTGCAAGTGACTCCAGCATCTTC
ACCN:	<u><a href="#">NM_138715</a></u> , <u><a href="#">NM_138715.1</a></u> , <u><a href="#">NM_138715.2</a></u> , <u><a href="#">BC063878</a></u> , <u><a href="#">BC063878.1</a></u> , <u><a href="#">NM_138715.3</a></u>
UniProt ID:	<u><a href="#">P21757</a></u>
Synonyms:	CD204; phSR1; phSR2; SCARA1; SR-A; SR-AI; SR-AII; SR-AIII; SRA
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)