

Product datasheet for **HP226780**

COMMD6 Human qPCR Primer Pair (NM_203497)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	170622
Forward Sequence:	ACCAGTGCCTCGTGTTCCTAG
Reverse Sequence:	GAAATGCTACCCTGCAAACAGGG
ACCN:	NM_203497 , NM_203497.1 , NM_203497.2 , NM_203497.3 , BC117391 , BC110837 , BC143912 , BC171758 , BM476717 , BX106312 , BX538116
UniProt ID:	Q7Z4G1
Synonyms:	Acrg
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)