

## Product datasheet for **HP226789**

### TCEAL4 Human qPCR Primer Pair (NM\_001006937)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	79921
Forward Sequence:	CAGCCACAAGACGAGAGAAAGC
Reverse Sequence:	TCCCTCACTCTCTGGCTTTCCT
ACCN:	<u><a href="#">NM_001006937</a></u> , <u><a href="#">NM_001006937.1</a></u> , <u><a href="#">NM_001006937.2</a></u> , <u><a href="#">BC012296</a></u> , <u><a href="#">BC012296.1</a></u> , <u><a href="#">BF359447</a></u> , <u><a href="#">BI084326</a></u> , <u><a href="#">BP375855</a></u> , <u><a href="#">BU536419</a></u> , <u><a href="#">BX104938</a></u>
UniProt ID:	<u><a href="#">Q96E15</a></u>
Synonyms:	NPD017; WEX7
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)