

Product datasheet for **HP226835**

FEZ2 Human qPCR Primer Pair (NM_005102)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	9637
Forward Sequence:	GTGATGCCTGTAGACTGGAAGTC
Reverse Sequence:	CACAGGAGACGATGATTGAGTGC
ACCN:	BC018032 , NM_005102 , NM_005102.1 , NM_005102.2 , BC018032.1 , BC110386 , BC110386.1 , BC037566 , BC047020 , BC054043 , BC064484 , BC089390 , BX537785
UniProt ID:	Q9UHY8
Synonyms:	HUM3CL
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)