

Product datasheet for **HP226864**

ORP1 (OSBPL1A) Human qPCR Primer Pair (NM_080597)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	114876
Forward Sequence:	AGAGCAGTCTCTGGTGAAAGGC
Reverse Sequence:	TCAAAGGAGCGTGCTGCTACTG
ACCN:	NM_080597 , NM_080597.1 , NM_080597.2 , NM_080597.3 , BC007004 , BC022293 , BC041563 , BC063420 , BC073854 , BC098384 , BC136452 , BC144117 , BX647893 , NM_080597.4
UniProt ID:	Q9BXW6
Synonyms:	ORP-1; ORP1; OSBPL1B
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)