

Product datasheet for **HP226869**

SULT1C2 Human qPCR Primer Pair (NM_176825)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	6819
Forward Sequence:	AGAAGTGCCAGCGAGCCATCAT
Reverse Sequence:	TTGAGACCAGCCTGGCAACAT
ACCN:	<u>BC005353</u> , <u>NM_176825</u> , <u>NM_176825.1</u> , <u>NM_176825.2</u> , <u>BC005353.1</u> , <u>BM849361</u> , <u>BT006951</u> , <u>NM_176825.3</u>
UniProt ID:	<u>O00338</u>
Synonyms:	humSULTC2; ST1C1; ST1C2; SULT1C1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)