

Product datasheet for **HP226921**

MRP6 (ABCC6) Human qPCR Primer Pair (NM_001171)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	368
Forward Sequence:	AGGCTTTCCTGCCCTTCTCCAT
Reverse Sequence:	CCAGAGGAACTTGAGTCTACGAC
ACCN:	<u>NM_001171</u> , <u>NM_001171.1</u> , <u>NM_001171.2</u> , <u>NM_001171.3</u> , <u>NM_001171.4</u> , <u>NM_001171.5</u> , <u>BC050733</u> , <u>BC131732</u>
UniProt ID:	<u>O95255</u>
Synonyms:	ABC34; ARA; EST349056; GACI2; MLP1; MOAT-E; MOATE; MRP6; PXE; PXE1; URG7
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)