

Product datasheet for **HP226975**

BTF (BCLAF1) Human qPCR Primer Pair (NM_014739)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	9774
Forward Sequence:	CCTAAACGAGCGGTTCACCTCG
Reverse Sequence:	GCTAAACGGGTATGCTTCCTCAG
ACCN:	<u>NM_014739</u> , <u>NM_014739.1</u> , <u>NM_014739.2</u> , <u>BC132780</u> , <u>BC016682</u> , <u>BC047687</u> , <u>BC047887</u> , <u>BC056894</u> , <u>BC063846</u> , <u>BC078160</u> , <u>BC144281</u> , <u>BC144282</u> , <u>BC144299</u> , <u>BX107174</u> , <u>BX537576</u> , <u>BX648431</u> , <u>BX648878</u> , <u>BX649160</u> , <u>NM_014739.3</u>
UniProt ID:	<u>Q9NYF8</u>
Synonyms:	bK211L9.1; BTF
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60 °C for 1 min; Stage 4: Melting curve: 95 °C for 15 sec, 60 °C for 15 sec, 95 °C for 15 sec.
Storage:	Store at -20 °C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)