

Product datasheet for **HP227000**

LIN54 Human qPCR Primer Pair (NM_194282)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	132660
Forward Sequence:	GTCCGACTTGTTACTGCCACATC
Reverse Sequence:	TTGGAAGTACATCCGCACTGG
ACCN:	<u>NM_194282</u> , <u>NM_194282.1</u> , <u>NM_194282.2</u> , <u>NM_194282.3</u> , <u>BC109277</u> , <u>BC109278</u> , <u>BX108765</u> , <u>BX537919</u> , <u>BX640657</u> , <u>BX640966</u>
UniProt ID:	<u>Q6MZP7</u>
Synonyms:	CXCDC1; JC8.6; MIP120; TCX1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)