

Product datasheet for **HP227004**

NEDD4 Human qPCR Primer Pair (NM_198400)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	4734
Forward Sequence:	CAGAAGAGGCAGCTTACAAGCC
Reverse Sequence:	CTCCCAACCTGGTGGTAATCC
ACCN:	<u>NM_198400</u> , <u>NM_198400.1</u> , <u>NM_198400.2</u> , <u>NM_198400.3</u> , <u>BC035707</u> , <u>BC047035</u> , <u>BC136605</u> , <u>BC144284</u> , <u>BC144285</u> , <u>BC152452</u> , <u>BC152562</u>
UniProt ID:	<u>P46934</u>
Synonyms:	NEDD4-1; RPF1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)