

Product datasheet for **HP227046**

UNC45B Human qPCR Primer Pair (NM_173167)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	146862
Forward Sequence:	CACCAATGGAGTGTCAGTCTC
Reverse Sequence:	GCAAACCTGCCTGAGACCGTAGT
ACCN:	<u>NM_173167</u> , <u>NM_173167.1</u> , <u>NM_173167.2</u> , <u>NM_173167.3</u> , <u>BC101062</u> , <u>BC101063</u> , <u>BX647642</u>
UniProt ID:	<u>Q8IWX7</u>
Synonyms:	CMYA4; CTRCT43; MFM11; SMUNC45; UNC-45B; UNC45
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)