

## Product datasheet for **HP227099**

### C2CD2 Human qPCR Primer Pair (NM\_199050)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	25966
Forward Sequence:	GCTGGTCTTTCATCAGTGTGCC
Reverse Sequence:	CAGAACCACTGATGGAGAGGCA
ACCN:	<u><a href="#">NM_199050</a></u> , <u><a href="#">NM_199050.1</a></u> , <u><a href="#">NM_199050.2</a></u> , <u><a href="#">BC020927</a></u> , <u><a href="#">BC035070</a></u> , <u><a href="#">BC062323</a></u> , <u><a href="#">NM_199050.3</a></u>
UniProt ID:	<u><a href="#">Q9Y426</a></u>
Synonyms:	C21orf25; C21orf258; TMEM24L
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)