

Product datasheet for **HP227120**

Syntabulin (SYBU) Human qPCR Primer Pair (NM_017786)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	55638
Forward Sequence:	AATCTGAGCGCCGACTCCATGA
Reverse Sequence:	CTTCCTGGCTTCTTTGAGTGCC
ACCN:	NM_017786 , NM_017786.1 , NM_017786.2 , NM_017786.3 , NM_017786.4 , NM_017786.5 , BC012082 , BG702715 , BU607867 , NM_017786.6
UniProt ID:	Q9NX95
Synonyms:	GOLSYN; OCSYN; SNPHL
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)