

Product datasheet for **HP227148**

ATPAF1 Human qPCR Primer Pair (NM_022745)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	64756
Forward Sequence:	AGGCATAGTGCTGATGACTGC
Reverse Sequence:	GGTCTCTTCCGATCAGTAGCG
ACCN:	<u>NM_022745</u> , <u>NM_022745.1</u> , <u>NM_022745.2</u> , <u>NM_022745.3</u> , <u>NM_022745.4</u> , <u>BC071731</u> , <u>BC071731.1</u> , <u>BC008498</u> , <u>BC018781</u> , <u>BX438342</u>
UniProt ID:	<u>Q5TC12</u>
Synonyms:	ATP11; ATP11p
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)