

Product datasheet for **HP227158**

MAGEA11 Human qPCR Primer Pair (NM_005366)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	4110
Forward Sequence:	CTCCAGTCAACAGAAAGAGCCC
Reverse Sequence:	CAGTGATCCTCTGGGTTCTCT
ACCN:	NM_005366 , NM_005366.1 , NM_005366.3 , NM_005366.4 , BC004479 , BC004479.2 , BQ212898 , BX337921 , NM_005366.5
UniProt ID:	P43364
Synonyms:	CT1.11; MAGE-11; MAGE11; MAGEA-11
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)