

Product datasheet for **HP227182**

ANKH Human qPCR Primer Pair (NM_054027)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	56172
Forward Sequence:	GTTACCTTCGTCTGCATGGCT
Reverse Sequence:	GGCAAAGTCCACTCCGATGATG
ACCN:	NM_054027 , NM_054027.1 , NM_054027.2 , NM_054027.3 , NM_054027.4 , BC014526 , BC014526.2 , BC009835 , BQ008340 , BU193060 , NM_054027.6
UniProt ID:	Q9HCJ1
Synonyms:	ANK; CCAL2; CMDJ; CPPDD; HANK; MANK; SLC62A1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)