

## Product datasheet for **HP227188**

### SUMF2 Human qPCR Primer Pair (NM\_015411)

#### Product data:

|                   |   |
|-------------------|---|
| Product Type:     | qPCR Primer Pairs   |
| Gene ID:          | 25870   |
| Forward Sequence: | GGAGACAAAGCTGAGGATGGCT  |
| Reverse Sequence: | TGGTACGGTGATGCTGTCCACT  |
| ACCN:             | <a href="#">BC008091</a> , <a href="#">NM_015411</a> , <a href="#">NM_015411.1</a> , <a href="#">NM_015411.2</a> , <a href="#">BC084539</a> , <a href="#">BC000224</a> , <a href="#">BC006159</a> , <a href="#">BC015600</a> , <a href="#">BC065222</a> , <a href="#">BC111092</a> , <a href="#">BM423367</a> , <a href="#">BP256067</a> , <a href="#">BP296294</a> |
| UniProt ID:       | <a href="#">Q8NBJ7</a>  |
| Synonyms:         | pFGE  |
| Components:       | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.  |
| Quality Control:  | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.                        |
| Storage:          | Store at -20°C.   |
| Stability:        | The primer mix is stable for one year from date of shipping.  |



[View online »](#)