

Product datasheet for **HP227197**

NSRP1 Human qPCR Primer Pair (NM_032141)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	84081
Forward Sequence:	AAAGAGGGCTGCTGCACTGGAA
Reverse Sequence:	GGCTTCACGAAAGCTGCATTTAG
ACCN:	NM_032141 , NM_032141.1 , NM_032141.2 , NM_032141.3 , BC040118 , BC016313 , BC046116 , BC062425 , BC105044 , BC105046 , BC143417 , BC143539 , NM_032141.4
UniProt ID:	Q9H0G5
Synonyms:	CCDC55; HSPC095; NSrp70
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)