

Product datasheet for **HP227240**

PTER Human qPCR Primer Pair (NM_030664)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	9317
Forward Sequence:	GGAGTTTGCTCAACTTGGCTGC
Reverse Sequence:	TCTTCACAGCCCTCTCCACCA
ACCN:	NM_030664 , NM_030664.2 , NM_030664.3 , NM_030664.4 , BC015092 , BC015092.1 , BC050411 , BC050411.1 , BU172463 , BX648556
UniProt ID:	Q96BW5
Synonyms:	HPHRP; RPR-1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)