

Product datasheet for **HP227275**

S100 Calcium Binding Protein A13 (S100A13) Human qPCR Primer Pair (NM_005979)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	6284
Forward Sequence:	GGAAGGATAGCCTCAGCGTCAA
Reverse Sequence:	CGAGTCCTGATTACATCCAAGC
ACCN:	<u>NM_005979</u> , <u>NM_005979.1</u> , <u>NM_005979.2</u> , <u>BC070291</u> , <u>BC068064</u> , <u>BC000632</u> , <u>BM762088</u> , <u>BM913375</u> , <u>BT006724</u> , <u>BX472819</u>
UniProt ID:	<u>Q99584</u>
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)