

Product datasheet for **HP227308**

SLC13A3 Human qPCR Primer Pair (NM_022829)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	64849
Forward Sequence:	CCATTGAGGAGTGGAACCTGCA
Reverse Sequence:	GTGTTGCTCAGCCACATGGACA
ACCN:	BC014931 , NM_022829 , NM_022829.1 , NM_022829.2 , NM_022829.3 , NM_022829.4 , NM_022829.5 , BC035966 , BC035966.1 , BM474324 , BP358260 , BQ946167 , NM_022829.6
UniProt ID:	Q8WWT9
Synonyms:	ARLIAK; NADC3; SDCT2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)