

Product datasheet for **HP227313**

DSN1 Human qPCR Primer Pair (NM_001145315)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	79980
Forward Sequence:	TACCCAGTGCTTCCAGAAGGTG
Reverse Sequence:	GTATGGCAGGTGGTTCTGTAG
ACCN:	<u>NM_001145315, NM_001145315.1, BC026011, BC035821, BC058899, NM_001145315.2</u>
UniProt ID:	<u>Q9H410</u>
Synonyms:	C20orf172; dj469A13.2; hKNL-3; KNL3; MIS13
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)