

## Product datasheet for **HP227343**

### DDX17 Human qPCR Primer Pair (NM\_006386)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	10521
Forward Sequence:	ACTGATGCAGCTTGTGGACCAC
Reverse Sequence:	AAGCCTTCGGTCACACTCATCC
ACCN:	<a href="#">NM_006386</a> , <a href="#">NM_006386.1</a> , <a href="#">NM_006386.2</a> , <a href="#">NM_006386.3</a> , <a href="#">NM_006386.4</a> , <a href="#">BC000595</a> , <a href="#">BC000595.2</a> , <a href="#">BC029553</a> , <a href="#">BC045747</a> , <a href="#">BI545691</a> , <a href="#">BP353248</a> , <a href="#">NM_006386.5</a>
UniProt ID:	<a href="#">Q92841</a>
Synonyms:	P72; RH70
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)