

Product datasheet for **HP227350**

TRAPPC5 Human qPCR Primer Pair (NM_001042462)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	126003
Forward Sequence:	TAGGCGCGTTGCTCTTCGCAA
Reverse Sequence:	GGAGATGTAGGTGTTGATGAGCG
ACCN:	<u>NM_001042462</u> , <u>NM_001042462.1</u> , <u>BC042161</u> , <u>BC042161.1</u> , <u>BC022836</u> , <u>BG528980</u> , <u>BI196928</u> , <u>BM841051</u> , <u>BU736689</u> , <u>NM_001042462.2</u>
UniProt ID:	<u>Q8IURO</u>
Synonyms:	TRS31
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)