

Product datasheet for **HP227411**

ADAM22 Human qPCR Primer Pair (NM_021723)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	53616
Forward Sequence:	CGGGTGCATAATGGGAGACACT
Reverse Sequence:	CCATTGCCACACTCAGGAGGAT
ACCN:	<u>NM_021723</u> , <u>NM_021723.1</u> , <u>NM_021723.2</u> , <u>NM_021723.3</u> , <u>NM_021723.4</u> , <u>BC125118</u> , <u>BC036029</u> , <u>BC062433</u> , <u>BC092465</u> , <u>BC125117</u> , <u>NM_021723.5</u>
UniProt ID:	<u>Q9P0K1</u>
Synonyms:	ADAM 22; DEE61; EIEE61; MDC2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)