

## Product datasheet for **HP227443**

### PSMA1 Human qPCR Primer Pair (NM\_148976)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	5682
Forward Sequence:	TGCAGAGCCATGTCCATTGGAG
Reverse Sequence:	GTAGTCAGGTCCTGTTCTGCAG
ACCN:	<a href="#">NM_148976</a> , <a href="#">NM_148976.1</a> , <a href="#">NM_148976.2</a> , <a href="#">BC005932</a> , <a href="#">BC005932.1</a> , <a href="#">BC002577</a> , <a href="#">BC008472</a> , <a href="#">BC009576</a> , <a href="#">BC015105</a> , <a href="#">BC015356</a> , <a href="#">BC022372</a> , <a href="#">BT006647</a>
UniProt ID:	<a href="#">P25786</a>
Synonyms:	HC2; HEL-S-275; NU; PROS30
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)